



## Endless Pilot Progressive Lenses

Fully personalized free-form progressive lens with a unique and innovative design that incorporates two zones for near vision.

For more information, contact

LAB LOGO  
HERE

**iot**

See the difference

**APPROVED PARTNER**

**###-###-####**

# Endless Pilot Progressive

The design architecture of **Endless Pilot Progressive** lenses is unique. In addition to a standard progressive configuration, it offers an extra segment for near vision at the top. It incorporates a lower vertical power progression with an upper addition segment. This creates areas for near vision at the top and bottom of the lens.



## Features & Benefits

- Fully personalized progressive lens
- 2 near vision areas: top & bottom
- Technology: IOT Digital Ray-Path 2
- Upper segment adapted to the wearer's visual needs
- Precise & comfortable near vision through the upper & lower areas of the lens
- Improved postural ergonomics avoiding unnecessary head movements
- Near elimination of peripheral blur
- Superior visual quality for viewing digital devices
- Excellent dynamic vision, easy transition between different viewing areas
- Comfortable & precise focus at all working distances
- MFHs: 16 & 18 mm